

Altitude XXL

Shell chair with castor wheels

VERMEIREN has designed **Altitude** with the **greatest comfort possible in mind for its user**. The Altitude chair in the XXL seat widths of 55 cm and 60 cm have a maximum load capacity of 150 kg. The upholstery of the shell ensures pressure relief and increases seating comfort, even during prolonged use. The inclination of the seat can be adjusted by an attendant and the angle-adjustable legrest can also be continuously adjusted by the user himself.




VERMEIREN
we care for you

Features

- 5 point belt
- Adjustable seat height
- Directional wheels
- Seat reclining up to 40°

Colours



- Removable safety belt
- Detachable, adjustable headrest



Option B52
Serum holder



Option:
378 table

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Altitude XXL** complies with the highest norms, laid down by certified test institutes in most European countries.

Technical data

Total width	800 mm
Total height	1280 - 1305 - 1330 mm
Total length	1150 mm
Seat width	550 / 600 mm
Seat height	570 - 595 - 620 mm
Seat angle	0° / 40°

Seat depth	510 mm
Backrest height	750 mm
Armrest height	230 mm
Weight	38 kg
Maximum user weight	150 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you