

EON

Inspired by the past,
turning to the future



ABOUT US

VERMEIREN, a manufacturer since 1957, is a leading name in the medical equipment industry. Over the years, the Vermeiren Group has specialised in producing wheelchairs, nursing beds, shell chairs, comfort chairs, and other products designed to support patients and users in maintaining their independence and mobility.

Vermeiren's European development and design teams use state-of-the-art precision tools, ensuring reliability and quality in every product they create. Today, Vermeiren is a global group with over 1,300 employees and subsidiaries across Europe, the UK and Asia, supplying products to more than 45 countries worldwide.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

Email : info@vermeiren.co.uk

Website : <https://www.vermeiren.co.uk>





 YouTube



SCAN AND
WATCH ME



DESIGNED WITH PRECISION

Our scooters are born from precision and pride, designed and built entirely in our own factories. Every element is thoughtfully crafted - no detail too small to matter. From the smooth feel of the steering column to the brightness of the LED lights, from the comfort of the adjustable chair to the quality of every nut and bolt, everything is scrutinised to

perfection. This isn't just about building a scooter; it's about creating a trusted companion for your daily journeys.

The result is a vehicle that blends function, safety, and style - something we are truly proud to put our name on, and you'll be proud to drive.



SBS* LED FRONT LIGHTS



LED INDICATORS



FOLDING BACKREST



ADJUSTABLE ARMRESTS



LED REAR LIGHTS



ANTI-TIPPING



DRIVE OR NEUTRAL MODE



CONTROLLER WITH BATTERY INDICATION



STEERING MUDGUARDS


*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.

COLOURS

OWN IT

Eon comes in a vibrant range of six colours, because your scooter should reflect you. Whether you want it to match to colour of your eyes, coordinate with your wardrobe, or simply choose your all-time favorite shade, colour is a powerful way to express your personality. A bold red might say "adventurous spirit," while a cool grey whispers "calm and confident."

From stylish neutrals to eye-catching brights, there's a perfect match for every mood and taste. After all, Eon isn't just transportation, it's part of your personal style, and the right color makes all the difference.



Arctic Green




Burgundy Red




Cool Grey



Matt Gold



Olive Green



Pearl Orange



Arctic Green

The Eon in Arctic Green reflects a calm, balanced personality with a love for nature and fresh beginnings - perfect for users who value harmony, optimism, and a refreshing, confident style.



Cool Grey

Why choose Cool Grey for your Eon? Because it speaks of modern confidence, understated elegance, and a calm, dependable nature, ideal for those who prefer style that whispers rather than shouts.



Matt Gold

Matt Gold reflects a personality that values refinement, confidence and timeless elegance. You appreciate life's finer details, and enjoy making a statement that feels both subtle and distinguished.






CUSTOMISE YOUR RIDE

Options transform a scooter from 'just transport' into your scooter. Mirrors, matching hubcaps and dust covers aren't just accessories—they're little touches that make every ride smoother, safer, and more enjoyable. A dust cover keeps it looking sharp.

Whether it's comfort, style, or convenience, customising your scooter lets you match it perfectly to your daily life. After all, why settle for standard when you can have a ride that's as unique, practical, and ready-for-anything as you are?

From left to right: crutch holder, set of mirrors, slick or studded tyres, coloured hubcap matching your scooter colour, shopping basket at the front, topcase at the back, storage protection cover.





Please contact your specialist dealer for information on compatibility.

TECHNIQUE

IN DEPTH

When it comes to scooters, the technical details aren't just numbers, they're the keys to freedom on four wheels. The width and length decide whether you glide gracefully through doorways or have to do a 2-point turn in the hallway. Weight tells you if it's a sturdy road companion or easier to nudge around.

Range means knowing if you can make it to the café and back without hunting for a charging plug. Even the top speed matters—fast enough to feel the breeze, slow enough to keep your hat on. These specifications make sure your ride is safe, comfy, and joyfully yours.

■ Total width	565 mm
■ Total length	1170 mm
■ Total height	1050 mm - 1100 mm
■ Folded height	900 mm
■ Seat height	460 mm - 520 mm
■ Seat width	470 mm
■ Seat depth	420 mm
■ Backrest height	480 mm
■ Footplate height	140 mm

■ Obstacle climbing	80 mm
■ Turning diameter	2700 mm
■ Batteries	38 Ah (C20), 12V DC, non-spillable AGM battery
■ Degree of protection	IPX4
■ Motors	320 W
■ Speed*	4 MPH or 6 MPH
■ Driving range*	19,9 or 18 Mile range
■ Weight	89 kg
■ Maximum user weight	140 kg

*The driving range is calculated according to European standards. *Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



VERMEIREN

we care for you



R.E.: **VERMEIREN GROUP**, Vermeirenp plein 1/15, 2920 Kalmthout - Belgium
info@vermeiren.com - www.vermeiren.com - BL Eon UK 2026-08 v1.0