

# Bretagne

## Relax chair

**Bretagne is a comfortable relaxing chair made by VERMEIREN. It can be placed in any room and makes for a fresh and homely atmosphere.**

This cosy chair has been developed with a luxurious touch. It is fully customizable to meet any taste and style, making it suitable for everyone.



**VERMEIREN**  
we care for you

## Comfort

Bretagne is available with or without wheels. The wheels have brakes to secure the chair into place. With this option, a push bar is mounted on the backrest as standard, so moving the Bretagne can be done in complete comfort and safety.

Bretagne can be obtained in several colours of wood and upholstery. These can be combined to your own liking. This way you can create the ideal chair for your interior.



## Products designed to suit you

**VERMEIREN** always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **Bretagne** complies with the highest norms, laid down by certified test institutes in most European countries.

### Upholstery



### Frame



## Technical data

■ Total width	620 mm
■ Total height	1180 mm (1260 mm*)
■ Total length	880 mm
■ Inclination backrest	13° - 45° mm
■ Seat width	500 mm
■ Height armrests	150 - 210 mm

■ Seat height	505 mm (555 mm*)
■ Height backrest	750 mm
■ Seat depth	480 mm
■ Weight	32,5 kg (34,75 kg*)
■ Maximum user weight	130 kg
■ Additional info	CE

\* With wheels. Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN GROUP

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer