

KUARZ

Inspired by the past,
turning to the future



ABOUT US

VERMEIREN, a manufacturer since 1957, is a leading name in the medical equipment industry. Over the years, the Vermeiren Group has specialised in producing wheelchairs, nursing beds, shell chairs, comfort chairs, and other products designed to support patients and users in maintaining their independence and mobility.

Vermeiren's European development and design teams use state-of-the-art precision tools, ensuring reliability and quality in every product they create. Today, Vermeiren is a global group with over 1,300 employees and subsidiaries across Europe, the UK and Asia, supplying products to more than 45 countries worldwide.

Email : info@vermeiren.com
Website : <https://www.vermeiren.com>





KUARZ



MOBILITY

WITH DIGNITY AND STYLE

A modern electric scooter is more than just a way to get from A to B - it's freedom, dignity, and style rolled into one. With its sleek, modern design, it's the kind of ride you're proud to be seen on, whether you're meeting friends for coffee, heading out for a bit of shopping, or enjoying a leisurely afternoon in the park. It's about moving through life with confidence,

knowing you have the independence to go wherever you want, whenever you want.

This isn't just mobility - it's a lifestyle upgrade, combining comfort, practicality, and a touch of flair that turns every trip into a little celebration of you.







DESIGNED WITH PRECISION

Our electric scooters are born from precision and pride, designed and built entirely in our own factories. Every element is thoughtfully crafted - no detail too small to matter. From the smooth feel of the steering column to the brightness of the LED lights, from the comfort of the adjustable chair to the quality of every nut and bolt, everything is scrutinised

to perfection. This isn't just about building a scooter; it's about creating a trusted companion for your daily journeys. The result is a vehicle that blends function, safety, and style - something we are truly proud to put our name on, and you'll be proud to drive.



SBS* LED FRONT LIGHTS



LED INDICATORS



COLLAPSIBLE CHAIR



PARKING BRAKE



LED REAR LIGHTS



ANTI-TIPPING



ROTATING AND ADJUSTABLE SEAT



LCD DISPLAY WITH ADAPTIVE BRIGHTNESS



PREMIUM SEAT FINISHING

*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.

COLORS

OWN IT

Kuarz comes in a vibrant range of six colors, because your scooter should reflect you. Whether you want it to match to color of your eyes, coordinate with your wardrobe, or simply choose your all-time favorite shade, color is a powerful way to express your personality. A bold red might say “adventurous spirit,” while a cool grey whispers “calm and confident.”

From stylish neutrals to eye-catching brights, there's a perfect match for every mood and taste. After all, Kuarz isn't just transportation, it's part of your personal style, and the right color makes all the difference.



Arctic Green

The Kuarz in Arctic Green reflects a calm, balanced personality with a love for nature and fresh beginnings - perfect for users who value harmony, optimism, and a refreshing, confident style.



Burgundy Red

Burgundy Red on the Kuarz scooter speaks of depth, warmth, and quiet strength. It's a choice for those who embrace life with passion, appreciate fine details, and exude an effortlessly classic charm.





Cool Grey

Why choose Cool Grey for your Kuarz? Because it speaks of modern confidence, understated elegance, and a calm, dependable nature, ideal for those who prefer style that whispers rather than shouts.



Olive Green

Olive Green as a color for your Kuarz? This embodies harmony, grounded confidence and a connection to nature, ideal for those who appreciate calm strength and a style that feels both fresh and enduring.



Mat Gold

Mat Gold reflects a personality that values refinement, confidence and timeless elegance. You appreciate life's finer details, and enjoy making a statement that feels both subtle and distinguished.



Pearl Orange

Pearl Orange suggests a personality that's energetic, creative and full of optimism. It reflects a zest for life, a love of standing out, and an ability to inspire warmth and positivity in others.



CUSTOMIZE YOUR RIDE

Options are what turn a scooter from “just transport” into your scooter. Mirrors, rearview cameras, and rain covers aren’t just accessories—they’re little touches that make every ride smoother, safer, and more enjoyable. A dust cover keeps it looking sharp, while a rollator holder means your mobility always comes with you.

Whether it’s comfort, style, or convenience, customizing your scooter lets you match it perfectly to your daily life. After all, why settle for standard when you can have a ride that’s as unique, practical, and ready-for-anything as you are?

From left to right: MOVO rain/sun protection, set of mirrors, slick or studded tyres, colored hubcap matching your scooter color, topcase at the back, shopping basket at the front, storage protection cover, rearview camera





TECHNIQUE

IN DEPTH

When it comes to scooters for seniors, the technical details aren't just numbers, they're the keys to freedom on four wheels. The width and length decide whether you glide gracefully through doorways or have to do a 2-point turn in the hallway. Weight tells you if it's a sturdy road companion or easier to nudge around.

Range means knowing if you can make it to the café and back without hunting for a charging plug. Even the top speed matters—fast enough to feel the breeze, slow enough to keep your hat on. These specifications make sure your ride is safe, comfy, and joyfully yours.

■ Total width	680 mm
■ Total length	1505 mm
■ Total height	1345 mm - 1395 mm
■ Folded height	870 mm
■ Seat height	480 mm - 530 mm
■ Seat width	525 mm
■ Seat depth	480 mm
■ Backrest height	550 mm
■ Footplate height	200 mm

■ Obstacle climbing	100 mm
■ Turning diameter	3400 mm
■ Batteries	80 Ah (C20), 12V DC, non-spillable AGM battery
■ Degree of protection	IPX4
■ Motors	750 W
■ Speed*	15 km/h
■ Driving range*	45 km
■ Weight	158 kg
■ Maximum user weight	150 kg

*The driving range is calculated according to European standards. *Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



KUARZ



SCAN AND
WATCH ME

PAST MEETS THE FUTURE

In our featured movie we strived to make a connection between the past and the future. Elegant lines and shapes are timeless, either in a classic mode of transportation that was built with love and care or in a new and trendy design for an electrical scooter. **Scan the QR code on top of this page and be prepared to discover the feel and looks of the Kuarz.**





VERMEIREN

we care for you



R.E. : **VERMEIREN GROUP**, Vermeirenplein 1/15, 2920 Kalmthout - Belgium
info@vermeiren.com - www.vermeiren.com - BL Kwarz World 2026-02 v1.2