

9062 - 9063

Commode chair

This classic commode chair with artificial leather upholstery delivers optimal comfort and hygiene. Sturdy fixed armrests and a strong backrest ensure the premium quality of this commode chair. The backrest is enveloping for a better seating comfort, the padded upholstery feels soft to the touch and the distance between the seat and the bucket is minimal to ensure better hygiene.



Model 9063
Height-adjustable legs



Model 9062





Enveloping backrest for a better seating comfort



9063 - Height adjustable legs




Soft padding for a softer body contact



Minimal distance between the seat and bucket

Frame colour(s) Upholstery colour(s)



C86



Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **9062 and 9063** comply with the highest norms, laid down by certified test institutes in most European countries.

Technical data

| | |
|-------------------|-------------------------------------|
| ■ Total width | 620 mm |
| ■ Total height | 9062 : 835 mm - 9063 : 755 - 885 mm |
| ■ Total length | 650 mm |
| ■ Seat width | 430 mm |
| ■ Height armrests | 215 mm |
| ■ Seat height | 9062 : 475 mm - 9063 : 470 - 560 mm |

| | |
|-----------------------|------------------------------|
| ■ Height backrest | 400 mm |
| ■ Seat depth | 450 mm |
| ■ Weight | 9062 : 8,35 kg - 9063 : 9 kg |
| ■ Maximum user weight | 120 kg |
| ■ Additional info | CE |

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you