



PLEGO

CARBON

LIGHT MOBILITY
WITHOUT COMPROMISE



MOBILITY EVERYDAY

Carbon frame, compact foldability, and everyday comfort. Perfect for active lifestyles, travel, and daily independence.

With its ultra-light carbon fibre frame, compact folding design, and premium comfort features, PLEGO CARBON is built for everyday independence. Whether you're at home, in the city, or travelling abroad, it offers the perfect balance of portability, performance, and style.

VISIT NOW
WWW.VERMEIREN.COM

MAKING SMOOTH
SHARP TURNS

FIT EASILY IN
YOUR CAR TRUNK

PERFECT FOR
YOUR TRAVELS

FREEDOM IN EVERY FOLD





CARBON

PREMIUM MATERIALS,
SUPERIOR PERFORMANCE

CRAFTED WITH A FULL CARBON FIBRE FRAME

STRENGTH, LIGHTNESS, AND ELEGANCE

EQUIPPED WITH MAGNESIUM ULTRA LIGHT
RIMS & FORKS

COMFORT FOR EVERY JOURNEY

PREMIUM BACKREST:
FOAM + 3D MESH



VISCOELASTIC
FOAM CUSHION



REMOVABLE
CUSHION



BACK
POCKET



OPTIONS

STANDING STRAIGHT
SOLUTION



ATTENDANT
CONTROLLER



EVERY JOURNEY, MADE COMFORTABLE

COMPACT FREEDOM,
MADE SIMPLE

FOLDS IN SECONDS

WE CARE FOR YOUR MOBILITY

One of the most compact powered wheelchairs in its class, PLEGO CARBON folds in seconds and fits easily into car trunk.

IDEAL FOR CITY LIFE,
TRAVEL, AND STORAGE



WE CARE FOR YOU

DESIGNED FOR CITY LIVING, TRAVEL,
AND EVERYDAY CONVENIENCE

TRAVEL MADE EASY

15,15 KG



Optional standing support

**TRAVEL WITH
PEACE OF MIND**

The Plego Carbon is equipped with a Lithium battery that complies with the international IATA* (International Air Transport Association) standards, giving you the freedom to fly with confidence.

* Please note: each airline has its own requirements regarding the transport of Lithium-ion batteries and retains the final decision on whether the wheelchair can be accepted on board. Before travelling, we recommend confirming with your airline or carrier that they allow battery-powered vehicles classified under UN3171. For further details, please refer to the user manual.

S M A R T PERFORMANCE

*DRIVING SETTINGS CAN BE ADJUSTED BY A QUALIFIED
PROFESSIONAL USING DEDICATED PROGRAMMING EQUIPMENT.

PLEGO CARBON delivers reliable performance for everyday life. The programmable* joystick lets you personalise acceleration, responsiveness, and control, giving you a truly tailored driving experience.



16 KM DRIVING RANGE WITH ONE BATTERY

PROGRAMMABLE JOYSTICK*

LINUX 250W MOTORS

12.8AH 25.9V LITHIUM BATTERY

WHY CHOOSE PLEGO CARBON?

Compared to similar folding carbon wheelchairs, PLEGO CARBON brings more than lightness. Its premium comfort, magnesium components, and compact design set a new standard. It's among the most compact powered chairs available, yet it also offers the comfort and adaptability that users need every day.



**LIGHTWEIGHT PORTABILITY, PREMIUM
COMFORT, AND SMART PERFORMANCE.**



More than lightness.
A smarter choice.

Category	Value
Width	590 mm
Length	915 mm
Seat width	450 mm
Seat depth	460 mm
Weight*	15,15 kg
Frame	Full carbon fibre (lightweight & strong)
Wheels & forks	Magnesium (light & durable)
Motors	Linux 250W
Max speed*	6 km/h
Drive range*	up to 16 km
Max. user weight	110 kg
Max kerb climb	15 mm
Turning radius	1030 mm
Battery	Lithium-ion 12.8Ah, 25.9V
Removable battery	Yes
Controller	Programmable joystick (personalised driving)
Backrest	Foam padding + breathable 3D mesh
Seat cushion	Visco foam with memory effect, removable
Folding	One of the most compact in its class
Transport	Easy to store and transport

* Weight excluding battery and seat cushion * The driving range is calculated according to European standards * Speed and driving range may deviate (ca. 15%).
Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°.
Subject to technical alterations.

CHECK FOR MORE
WWW.VERMEIREN.COM

PLEGO

CARBON



**LIGHTWEIGHT. COMPACT.
PREMIUM COMFORT.**



VERMEIREN GROUP

Vermeirenplein 1-15

B-2920 Kalmthout

Belgium

www.vermeiren.com

info@vermeiren.com