



# SYNC

Inspired by the past,  
turning to the future



# ABOUT US

**VERMEIREN**, a manufacturer since 1957, is a leading name in the medical equipment industry. Over the years, the Vermeiren Group has specialised in producing wheelchairs, nursing beds, shell chairs, comfort chairs, and other products designed to support patients and users in maintaining their independence and mobility.

Vermeiren's European development and design teams use state-of-the-art precision tools, ensuring reliability and quality in every product they create. Today, Vermeiren is a global group with over 1,300 employees and subsidiaries across Europe, the UK and Asia, supplying products to more than 45 countries worldwide.

## **VERMEIREN** GROUP

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Email : [info@vermeiren.com](mailto:info@vermeiren.com)

Website : <https://www.vermeiren.com>





**PURA BRASA**

**PESCADOS** **CARNES**

**PASTA FRESCA**

**PURA SURGIDOS** **CARNES A PIEDRA**



YouTube



SCAN AND WATCH ME



# DESIGNED WITH PRECISION

Our scooters are born from precision and pride, designed and built entirely in our own factories. Every element is thoughtfully crafted - no detail too small to matter. From the smooth feel of the steering column to the brightness of the LED lights, from the comfort of the adjustable chair to the quality of every nut and bolt, everything is scrutinised to

perfection. This isn't just about building a scooter; it's about creating a trusted companion for your daily journeys.

The result is a vehicle that blends function, safety, and style - something we are truly proud to put our name on, and you'll be proud to drive.



SBS\* LED FRONT LIGHTS



LED INDICATORS



COMPACT STORAGE



ROTATING AND ADJUSTABLE SEAT



LED REAR LIGHTS



ANTI-TIPPING



DRIVE OR NEUTRAL MODE



CONTROLLER WITH BATTERY INDICATION



FOLDING BACKREST

\*SBS: See and Be Seen - Advanced LED light technology for increased safety when driving in low light conditions.

# COLOURS

## OWN IT

---

Sync comes in a vibrant range of six colours, because your scooter should reflect you. Whether you want it to match to colour of your eyes, coordinate with your wardrobe, or simply choose your all-time favorite shade, colour is a powerful way to express your personality. A bold red might say “adventurous spirit,” while a cool grey whispers “calm and confident.”


From stylish neutrals to eye-catching brights, there's a perfect match for every mood and taste. After all, Sync isn't just transportation, it's part of your personal style, and the right colour makes all the difference.



Arctic Green



Burgundy Red




Cool Grey



Matt Gold



Olive Green



Pearl Orange



### Burgundy Red

Burgundy Red on the Sync scooter speaks of depth, warmth, and quiet strength. It's a choice for those who embrace life with passion, appreciate fine details, and exude an effortlessly classic charm.



### Cool Grey

Why choose Cool Grey for your Sync? Because it speaks of modern confidence, understated elegance, and a calm, dependable nature, ideal for those who prefer style that whispers rather than shouts.



### Olive Green

Olive Green as a colour for your Sync? This embodies harmony, grounded confidence and a connection to nature, ideal for those who appreciate calm strength and a style that feels both fresh and enduring.





# CUSTOMISE YOUR RIDE

---

Options transform a scooter from 'just transport' into your scooter. Mirrors, matching hubcaps and rain covers aren't just accessories—they're little touches that make every ride smoother, safer, and more enjoyable. A dust cover keeps it looking sharp, while a rollator holder means your mobility always comes with you.

Whether it's comfort, style, or convenience, customising your scooter lets you match it perfectly to your daily life. After all, why settle for standard when you can have a ride that's as unique, practical, and ready-for-anything as you are?

**From left to right:** MOVO rain/sun protection, set of mirrors, slick or studded tyres, coloured hubcap matching your scooter colour, topcase at the back, shopping basket at the front, storage protection cover, Quadri-Light rollator holder, 286i rollator holder, crutch holder





# TECHNIQUE

## IN DEPTH

When it comes to scooters, the technical details aren't just numbers, they're the keys to freedom on four wheels. The width and length decide whether you glide gracefully through doorways or have to do a 2-point turn in the hallway. Weight tells you if it's a sturdy road companion or easier to nudge around.

Range means knowing if you can make it to the café and back without hunting for a charging plug. Even the top speed matters—fast enough to feel the breeze, slow enough to keep your hat on. These specifications make sure your ride is safe, comfy, and joyfully yours.

■ Total width	630 mm
■ Total length	1235 mm
■ Total height	1275 mm - 1325 mm
■ Folded height	830 mm
■ Seat height	420 mm - 470 mm
■ Seat width	470 mm
■ Seat depth	495 mm
■ Backrest height	530 mm
■ Footplate height	185 mm

■ Obstacle climbing	110 mm
■ Turning diameter	3530 mm
■ Batteries	38 Ah or 50 Ah (C20), 12V DC, non-spillable AGM battery
■ Degree of protection	IPX4
■ Motors	450 W / 500 W / 650 W
■ Speed*	6 km/h / 10 km/h / 12 km/h / 15 km/h
■ Driving range*	28 km - 40 km
■ Weight	109 kg
■ Maximum user weight	140 kg

\*The driving range is calculated according to European standards. \*Speed and driving range may deviate (ca. 15%). Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



# VERMEIREN

we care for you



R.E. : **VERMEIREN** group, Vermeirenplein 1/15 - 2920 Kalmthout - Belgium  
info@vermeiren.com - www.vermeiren.com - BL Sync World 2026-05 v1.1