




Version: August 2024

All rights reserved, including translation. No part of this manual may be reproduced in any form what so ever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated or distributed by using electronic systems.

**VERMEIREN** GROUP



## Kids

2202 Aqua - Safari

2215 Happy - Sporty

6

6

10

## Adolescents

2217 Freedom

14

14

## Adults

2219 Lagoon

2219 Vintage

E-version

16

16

18


20

## Options

22

## How do I choose my tricycle?

The choice of your tricycle is determined by the total length of the inner leg (A).  
(total distance, stretched leg, from the folds of the groin to the heel)





2202



Inside leg  
from 30 to 46 cm

3 years  
and  
above

## AQUA - unisex model

Decorative figures and bright, shiny colours form the main attraction for young children. The bell is shaped as a dolphin and the rear tray is ideal to transport toys around. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



## Range of options available at page 22.

								
total length	total width	saddle height	height saddle-pedal	step-in height	height steering	wheel size	total weight	users weight
1110	648	500-640	390-520	118	726-838	12" x 2.00	14	40

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.



002

# SAFARI

2202

Inside leg  
from 30 to 46 cm

3 years  
and  
above

## SAFARI - unisex model

Decorative safari figures and bright colours make this tricycle eye-catching to young children. The bell is giraffe shaped and a rear tray is ideal to transport toys. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



## Range of options available at page 22.

									
total length	total width	saddle height	height saddle-pedal	step-in height	height steering	wheel size	total weight	kg	users weight
1110	648	500-640	390-520	118	726-838	12" x 2.00	14	kg	40

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.



# Happy

## HAPPY - 100% girls

Cute pink colours in combination with a clear white and a rear tray, ideal to transport toys. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



8 years and above



V-brakes

2215

Inside leg  
from 50 to 70 cm

### Range of options available at page 22.



total  
length

1294



total  
width

648



saddle  
height

610-760



height  
saddle-pedal

48-650



step-in  
height

133



height  
steering

827-940



wheel  
size

16" x 1,75



kg

total  
weight

16,5



kg

users  
weight

60

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.

2215



# Sporty

2215

Inside leg  
from 50 to 70 cm

8 years  
and  
above

## SPORTY - Boys will be boys

Deep dark black topped with fiery red flames and a rear tray, ideal to transport things. The steering wheel (different steerings available) is adjustable in height and angle to make this tricycle a perfect fit to the needs of your child.

- Handbrake with parking lock.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- Fixed pinion for continuous movement (forward or backwards).
- Different options available.



## Range of options available at page 22.



total  
length

1294



total  
width

648



saddle  
height

610-760



height  
saddle-pedal

48-650



step-in  
height

133



height  
steering

827-940



wheel  
size

16" x 1,75



total  
weight

16,5



users  
weight

60

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5") are possible.

2217

# FREEDOM

## WHITE

The Freedom is a sturdy tricycle ideal for demanding adolescents available in black and white. A perfect tricycle for strolling through countryside.

Range of options available at page 22. Also available as e-Bike at page 20.



total length

1599



total width

670



saddle height

720-860



height saddle-pedal

590-780



step-in height

174



height steering

986-1174



wheel size

20" x 1,75



total weight

25,5



users weight

110

# FREEDOM BLACK

- Different options available.
- The steering wheel (different steerings available) is adjustable in height and angle.
- Saddle adjustable in height.
- Lowered frame for easy mounting/dismounting.
- V brakes at the front with parking lock. Optional disc brake for the rear wheels available.
- Big metal transport tray at the rear.



2217

Inside leg  
from 68 to 88 cm

14 years  
and  
above



V-brakes

**LAGOON - Adult tricycle**

Available in a classic bordeaux colour or a more neutral grey tone. The perfect tricycle for long rides, offering stability and comfort.

**Range of options available at page 22.**

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



# L A G O O N

- The steering wheel (different steerings available) is adjustable in height and angle.
- V brakes at the front or an optional disc brake for the rear wheels.
- Lowered frame for easy mounting/dismounting.
- Saddle adjustable in height.
- Big metal transport tray at the rear.



Inside leg  
from 75 to 108 cm

Adults



V-brakes



total  
length  
1808



total  
width  
768



saddle  
height  
810-900



height  
saddle-pedal  
640-880



step-in  
height  
194



height  
steering  
1121-1309  
1013-1206 (Vintage)



wheel  
size  
24" x 1<sup>3/8</sup>



kg

total  
weight  
27,8



kg

users  
weight  
125

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

**VINTAGE - Adult tricycle**

A classic and stylish retro will get the peoples attention. Brown tyres and a large optional vintage basket completes the picture.

**Range of options available at page 22.**

Also available in **E-version**, with an integrated electric motor which can be used for propulsion.



# Vintage

2219

Inside leg  
from 75 to 108 cm

Adults

- The steering wheel (different steerings available) is adjustable in height and angle.
- V brakes at the front with parking lock or an optional disc brake for the rear wheel.
- Lowered frame for easy mounting/dismounting.
- Big metal transport tray at the rear. Optional vintage basket.
- Saddle adjustable in height.



total  
length  
1808



total  
width  
768



saddle  
height  
810-900



height  
saddle-pedal  
640-880



step-in  
height  
194



height  
steering  
1121-1309  
1013-1206 (Vintage)



wheel  
size  
24" x 1<sup>3/8</sup>



total  
weight  
27,8



users  
weight  
125

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

2217E



Manual controller



Battery pack

# E-version

Electrical model with propulsion support.

The E-version for adults is available for the **LAGOON** and the **VINTAGE** version. The electric tricycle will allow you to go further, easier and faster.

The operating system is very easy to use. A simple **ON/OFF** switch at the battery and a controller at the steering wheel to get started. It is also equipped with a **BOOSTER** button to give a quick acceleration to the tricycle (for example when you need to start uphill).

The battery pack, located at the back of the tricycle is light and removable, allowing you to easily recharge it at home.

The electrical assistance will power the front wheel of the tricycle.



E-Version

Inside leg  
from 75 to 108 cm

Adults



17 km/h



36 V / 9 Ah



250 W

Autonomy : depends on the usage, state of the road, body features and selected propulsion.

# COMPATIBILITY TABLE

		2202	2215	2217	2219	2219 vintage
<b>S1</b>	Classic steering	•	•	•	•	-
<b>S2</b>	Flat steering	•	•	•	•	-
<b>S3</b>	City steering	-	-	•	•	standard
<b>S4</b>	Rectangular steering	•	•	•	-	-
<b>S5</b>	Adjustable steering bracket	-	-	•	•	•
<b>D1</b>	Steering limiter	•	•	•	•	•
<b>06</b>	Fixed wheel	•	•	•	•	•
<b>01 + B1</b>	Free wheel	-	-	•	•	•
<b>B1</b>	Disc brake for the rear heel	-	-	•	•	•
<b>02</b>	Torpedo brake	•	•	•	•	•
<b>03</b>	3 gears	-	•	•	•	•
<b>07</b>	7 gears	-	-	•	•	•
<b>23</b>	Differential	-	-	-	•	•
<b>F1</b>	Pedals with strap	•	•	•	•	•
<b>F2</b>	Pedals with toe clip & strap	•	•	•	•	•
<b>F3</b>	Foot support with Velcro straps	• (S)	• (S)	• (M)	• (M)	• (M)
<b>F4</b>	Foot support with Velcro straps & ankle supp.	•	•	•	•	•
<b>F5</b>	Pedals with adjustable length	-	•	•	•	•
<b>E1</b>	Pushing handle	• (S)	• (S)	• (M)	• (M)	• (M)

		2202	2215	2217	2219	2219 vintage
<b>E2</b>	Lighting system	-	•	•	•	•
<b>E3</b>	Odometer	-	•	•	•	•
<b>E4</b>	Vintage basket	-	-	-	-	•
<b>E5</b>	Lockable topcase	-	-	•	•	•
<b>E6</b>	Crutch holder	-	-	•	•	•
<b>E7</b>	Drink holder	•	•	•	•	•
<b>E8</b>	Mirror	•	•	•	•	•
<b>E9</b>	Hand straps	•	•	•	•	•
<b>E10</b>	Hands holder with straps	• (S)	• (S)	• (M)	• (M)	• (M)
<b>E11</b>	Closed hands holder	•	•	•	•	•
<b>E12</b>	Comfort saddle	-	-	•	•	•
<b>E13-D</b>	Home trainer	-	•	•	•	•
<b>E14</b>	Front basket	-	-	•	•	•
<b>BAT 9</b>	Battery pack E-bike	-	-	-	•	•
<b>P1</b>	Pelvic support	• (S)	• (S)	• (M)	• (M)	• (M)
<b>P2</b>	Back support	• (S)	• (S)	• (M)	• (M)	• (M)
<b>P3</b>	Pelvic + back support (P1 + P2)	• (S)	• (S)	• (M)	• (M)	• (M)
<b>P4</b>	Adjustable saddle	• (S)	• (S)	• (M)	• (M)	• (M)
<b>P5</b>	Saddle tube with shock absorber	-	-	-	-	•
<b>E-Tricycle</b>	With electric pedal support	-	-	•	•	•

(S) : SMALL - (M) : MEDIUM

# CONFIGURATION

## Steering



**S1** - Classic steering



**S2** - Flat steering



**S3** - City steering



**S4** - Rectangular steering



**S5** - Adjustable steering bracket (-10° to +50°)



**D1** - Steering limiter

## Moving



**O6** - Fixed wheel



**O1 + B1** - Free wheel



**B1** - Disc brake for the rear wheel



**O2** - Torpedo brake



**O3** - 3 gears  
**O7** - 7 gears



**E13** - Hometrainer

## Support



**P1** - Pelvic support



**P2** - Back support



**P3** - Pelvic and back support



**E9** - Hand straps



**E10** - Hands holder with straps



**E11** - Closed hands holder

# A N D O P T I O N S

## Miscellaneous



**BAT 9** - Extra battery  
E-bike 9A



**P4** - Adjustable saddle  
(140 mm)



**P5** - Saddle tube with  
shock absorber



**E12** - Comfort saddle  
PU + plastic



**E1** - Pushing handle  
(335 mm or 605 mm)



**E14** - Front basket

## Pedals



**F1** - Pedals with strap



**F2** - Pedals with toes  
clip & strap



**F3** - Foot support  
with Velcro straps



**F4** - Foot support  
with Velcro straps  
& ankle support



**F5** - Pedals with  
adjustable length  
(22 mm)



**E2** - Lighting system  
(batteries included)

## Accessories



**E3** - Odometer



**E4** - Vintage Basket



**E5** - Lockable topcase  
volume: 30 liter



**E6** - Crutch holder



**E7** - Drink holder  
(cannister not included)




**E8** - Rear mirror

**VERMEIREN GROUP**

Vermeirenplein 1-15  
B-2920 Kalmthout

Tel: +32 (0)3 620 20 20  
Fax: +32 (0)3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



**VERMEIREN**  
we care for you