


* - For Ecl. X4 Kids use small pads:
XC00686 (3992455) /
Alternative: XP03265/2996861

- Only left side depicted
- Quantities for one side

No.	Plan	Article	Qty
1	E001038	1913465 (R) 1913464 (L)	1 1
2	P00046914	2997543	1
3	M6 DIN125A	1001564	2
4	M5 DIN985	1001679	4
5	M5x20 DIN912	1003907	4
6	M6x12 DIN986	1008492	2
7	M6x55 DIN912	1908726	2

Method PL		
Procedure	Description	Tool
Drill	Drilling 2 holes $\phi 6$	M500052

** - Drilling instruction:



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

BI5 ARMREST ECLIPS X4 KIDS				---
Scale: 10:33	Drafted	Verified	Confirmed	Plan No. E00116355/A.1
Date	2.10.2025	02/10/2025	02/10/2025	
Initials	AZ	AZ	AS	