



Bill of Material			
No.	Plan	Article	Quantity
1	M8 DIN985	1996213	1
2	X002133	1992424	4
3	X003692	1998763	4
4	X004452	1904058	4
5	X004606	1905091 (39/44/48)	2
	X004607	1905092 (42/46/50)	2
6	X004620	1905090	2
7	P00023336	2639354 (39)	1
	P00023305	2639554 (42)	1
	P00023279	2638154 (44)	1
	P00023307	2639754 (46)	1
	P00023309	2639954 (48)	1
	P00023311	2640154 (50)	1
8	P00023337	2639454 (39)	1
	P00023306	2639654 (42)	1
	P00023280	2638254 (44)	1
	P00023308	2639854 (46)	1
	P00023310	2640054 (48)	1
	P00023312	2640254 (50)	1
9	3,9x9,5 DIN7981	2992349	4
10	M8 DIN126	1996214	2
11	M4 DIN985	2994901	4
12	M6 DIN985	2994063	6
13	M6x40 ISO7380	2994066	4
14	M6x50 DIN912	2995258	2
15	M8x105 DIN912	2997555	1
16	M4x16	2995260	4
17	P000064	1905711	1
18	P000064	1916218	1

CROSS / V500 ACTIVE				-
				Plan No. E00050163/C1.1
Scale: 1:5	Drafted	Verified	Confirmed	
Date	1/09/2021	19/09/2024	19/09/2024	
Initials	PM	AZ	AS	